



Respimatic 100



RESPIMATIC 100

(Patent Pending)

Is it right for you?

Only non-invasive support?

Support from Initiation to Weaning?

No compressed air or piped Oxygen?

Connect to O₂ Cylinder or Concentrator?

Full range of Respiration parameters?

Breath Synchronization for Patient Comfort?

Remote monitoring capability?

Handle harsh-uncontrolled Environment?

Easy-to-use System?


Budget Friendly?

Setting the Context

The Motivation


The Problem Statement

The Motivation Serve the Bottom of the Pyramid



The Problem

2 BIG Challenges




Lack of Affordability leading
to Scarcity



Lack of Skilled Practitioners

Observations on Ventilator Evolution



Iron Lung
Age




Pneumatic
Age



μController
Age



Smart “E”-
Age



Most-used Ventilation Modes have not changed

- Volume and Pressure Control
- Control BPM, I/E, VT and PS
- Monitor pressures and flow
- Safety Alarm systems

Diminishing Returns from what has evolved ...

- Exotic Ventilation modes
- Multitude of Sensors
- Fancy Touch-screen LCD Displays

Respimatic – Respiratory Technology Revisited

Remote Diagnosis and Monitoring (Telemedicine)

Rural, remote geographies do not have skilled specialists
BUT Communication Technologies have come of age

Leverage Time- proven Technologies

Simplify design and feature list
Simplify HMI
Use Off-the-shelf proven components

Maintenance Breaths in case of unexpected errors

System cannot stop delivering breaths.
Implement appropriate **Backup** mechanism for each scenario.



Respimatic 100 System Details

System Components

Technical Details

Our Solution **RESPIMATIC 100**

4 Commonly Used
Ventilation Modes
CMV, ACV, SIMV, PSV

Respiration Rate, Tidal
Volume, PEEP, Support
Pressure Controls

Volume Controlled and
Pressure Supported
Breaths

Mandatory &
Spontaneous Breaths
with Full Breath
Synchronization

Complete set of WEB Apps

Remote Dashboard
Remote Recorder
Remote Analyzer
FiO₂ Calculator

Low-speed Wi-Fi sufficient
Phone Hot-spot also sufficient

Uses secure HTTPS protocol

RESPIMATIC 100 WEB APPS

Calculate FiO₂ Settings

Select a Known System

BANGALORE

Add a Known System

Manage Known Systems

Launch Dashboard

Launch Analyzer

Delete ALL Cookies



Right icons created by Freepik - Flaticon

Remote WEB Dashboard Snapshot View


Anyone, anywhere in the world can monitor any patient via the WEB

Multiple specialists can monitor the same patient

One specialist can monitor multiple patients

5 Dashboard views

- Snapshots
- Charts
- Statistics
- Breath Shapes
- Alerts



<https://www.respimatic.com>

Remote WEB Dashboard – Charts View



Remote WEB Dashboard – Statistics View

Shape Session [20-12-2022 10:53:21]

↻ ↻ ✓ ✕

Parameter	Units	Min	Max	Avg
Peak Pressure	cmH20	1	35	28.1
Plateau Pressure	cmH20	17	33	27.1
PEEP Pressure	cmH20	5	7	6.0
Tidal Volume Delivered	ml	150	412	359.4
Minute Volume Delivered	litres/min	5.6	5.9	5.7
Mandatory BPM	bpm	15	16	15.2
Spontaneous BPM	bpm	0	1	0.7
FIO2	%	50	50	50.0
Instantaneous Static Compliance	ml/cmH20	14	30	18.4
Instantaneous Dynamic Compliance	ml/cmH20	14	22	16.5
System Temperature	degC	29	30	29.2


Parameter	Units	Values
Ventilation Mode	mode	ACV,CMV
Tidal Volume	ml	400,300
Respiration Rate	bpm	15,20
I:E Ratio	ratio	1:3,1:2
PEEP Pressure	cmH20	6
Maximum Pressure	cmH20	50
Support Pressure	cmH20	25
Support Pressure Termination	%flow,secs	F20%
FIO2	%	50

Information	Value
Number of Breaths	73
Number of Mandatory Breaths	56
Number of Spontaneous Breaths	17
Number of Maintenance Breaths	0
Number of Missing Breaths (Comms Failure)	0
Number of Notifications	0
Number of Warnings	2
Number of Errors	17

MODE	VT	RR	I:E	PEEP	PPMAX	PS	TPS	FIO2	# of BREATHS	Before BREATH#
ACV	400	15	1:3	6	50	25	F20%	50	17	0
CMV	300	15	1:2	6	50	25	F20%	50	16	17
ACV	400	15	1:3	6	50	25	F20%	50	21	33
ACV	400	20	1:3	6	50	25	F20%	50	19	54

Respimatic 100 - Preliminary and Confidential

Remote WEB Dashboard – Breath Shapes View



Remote WEB Dashboard – Alerts View

The screenshot displays a dashboard interface with three main sections: Session Errors, Session Warnings, and Session Information. Each section contains a list of alerts with their respective dates and times.

Session Errors:

- ERROR #1 DateTime: [09-12-2022]06:55:16
Leakage in
Breathing Circuit
Switching to
Maintenance Breaths
- ERROR #2 DateTime: [09-12-2022]06:55:21
[ERROR] state
Press PAUSE to show
the System state
leading to ERROR
- ERROR #3 DateTime: [09-12-2022]06:55:25
Leakage in
Breathing Circuit
Maintenance mode
Deliver safe breaths
- ERROR #4 DateTime: [09-12-2022]06:55:29
[ERROR] state

Session Warnings:

- WARNING #1 DateTime: [09-12-2022]06:43:16
PEEP delta measured
up to -0.8 cm H2O
Adjust valve/setting
YES -> Commit
- WARNING #2 DateTime: [09-12-2022]06:43:19
PEEP delta measured
up to -1.0 cm H2O
Adjust valve/setting
YES -> Commit
- WARNING #3 DateTime: [09-12-2022]06:43:21
PEEP delta measured
up to -0.9 cm H2O
Adjust valve/setting
YES -> Commit
- WARNING #4 DateTime: [09-12-2022]06:55:14
PEEP delta measured

Session Information:

- INFO #1 DateTime: [09-12-2022]06:45:23
1 Breath(s) missed
Info not received by
Dashboard due to
Internet packet loss
- INFO #2 DateTime: [09-12-2022]06:46:26
1 Breath(s) missed
Info not received by
Dashboard due to
Internet packet loss
- INFO #3 DateTime: [09-12-2022]06:47:17
1 Breath(s) missed
Info not received by
Dashboard due to
Internet packet loss
- INFO #4 DateTime: [09-12-2022]06:47:28
1 Breath(s) missed

At the top of the dashboard, there is a header bar with the number "1" on the left, a "RESET" button with a hand icon, a "Set Breath Range Window for Display" checkbox with a checked status and a checkmark icon, and a close button with an "X" icon on the right. The total count "204" is displayed at the top right.

Remote WEB Analyzer

Any patient Session can be recorded locally or remotely.

The recorded Session can then be analyzed off-line using the Analyzer.

Demo Session [09-12-2022 06:39:24]

RESPIMATIC-100 Session Databases		
Session Name	Created	Actions
Demo Session	09-12-2022 06:39:24	<input checked="" type="checkbox"/> <input type="button" value="Edit"/> <input type="button" value="Delete"/>
New Session	29-11-2022 09:08:55	<input checked="" type="checkbox"/> <input type="button" value="Edit"/> <input type="button" value="Delete"/>
Demo	29-11-2022 07:51:30	<input checked="" type="checkbox"/> <input type="button" value="Edit"/> <input type="button" value="Delete"/>

SYSTEM UID


UID_28615E07D6013C4A
(BANGALORE)

RESPIMATIC 100 WEB ANALYZER



<https://www.respimatic.com>

WEB Apps Architecture



The Human-Machine Interface – Front Panel

Simple Tactile buttons
No delicate touch screen etc.

Easy to read 7-seg LED Parameters Display

Parameter selection using simple arrow buttons

4-line LCD Display for displaying Messages and Menus

Peak, Plateau, PEEP pressures displayed after each breath

Also shows Delivered Volumes, Lung Compliance, Breath types etc.




System Architecture

Clearly separated
Frontend, Backend and
Breathing system

IP is in the Frontend
design and Algorithms

Backend and Breathing
System can have
multiple avatars



Respimatic 100

Under the hood

Production Cost Rs. 50,000

Potentially lower with bulk negotiations

Compact, Lightweight, and Robust

Usable in harsh environments

Intuitive HMI - Simple to operate

Front Control Panel & Remote Monitoring





Simple electronics
with time proven
COTS components

Rugged mechanical
system – various
alternatives possible

Essential Parameter
monitoring and
carefully vetted
alarm conditions

Robust and easy to
clean enclosure
surfaces

Suitable for mass
production



so, $(P_{G1} - P_{G2})$ can be replaced by $E * (P_{G1} - P_{G2})$ for an appropriate transforms to the one below.

$$Q = K * \sqrt{E * (P_{G1} - P_{G2})}$$

$$\text{where } K = \frac{C}{\sqrt{\rho \text{atmosphere} * \sqrt{E}}}$$

$$\text{or } K = \frac{f(Re)}{\sqrt{\rho \text{atmosphere} * \sqrt{E}}} \text{ where } Re \text{ is the Reynolds number}$$

Rewriting for every sample interval time t , the equation is as below.

$$Q(t) = K * \sqrt{E * (P_{G1}(t) - P_{G2}(t))}$$

Breathing Circuit

One proprietary, **patent-pending** Dual Pressure line connector with Orifice plate

Off-the-shelf single-limb Breathing Circuit with NRBV

BVM or Ambu Bag with Reservoir

Pressure sensors, PEEP valve

HME Filter


Humidifier

Oxygen Source



Off-the-Shelf Single limbed Circuit with
NRBM

Patent-pending Dual
Pressure line Connector



Ventilation Modes

The 4 most frequently used

<i>Continuous Mandatory Ventilation (CMV)</i>	<i>Synchronized Assist Control Ventilation (Sync ACV)</i>	<i>Synchronized Intermittent Mandatory Ventilation (SIMV)</i>	<i>Pressure Support Ventilation (PSV)</i>
Volume Controlled Mandatory Breaths Ignore spontaneous breaths	Volume Controlled Mandatory Breaths Volume controlled breaths in response to spontaneous breaths Breath Synchronization	Volume Controlled Mandatory Breaths Pressure supported breaths in response to spontaneous breaths Breath Synchronization	Pressure supported breaths in response to spontaneous breaths Careful Monitoring of Minute Volume No Mandatory breaths except when in dire need

Volume Controlled Breaths

(All modes)

Tidal Volume (ml)

200 to 600 ml
increments of 50 ml

Respiratory Rate (bpm)

10 to 30 bpm
increments of 1

Inspiration/Expiration Ratio (I:E)

1:1, 1:2, 1:3

PEEP (cmH₂O)

4 to 15 cmH₂O
increments of 1 cmH₂O

Max Pressure (cmH₂O)

20 to 60 cmH₂O
increments of 5 cmH₂O

FiO₂ Support

Externally Table Driven
21% to 100%

Pressure Supported Breaths

(SIMV & PSV modes)

Support Pressure (PS)

10 cm H₂O to 40 cm H₂O in increments of 5 cm H₂O

Support Pressure Termination (TPS)

Flow-dependent

Terminate when flow falls below 10%, 20% or 30%
of peak flow

Time dependent

Terminate after 1.0 secs to 3.0 secs in increments of
0.5 secs

Both ACV and SIMV modes

- A must for patient comfort
- Synchronize Mandatory breaths with Spontaneous breaths
- Prevent breath stacking

Breath Synchronization

FiO_2 Settings


Without the Reservoir bag, FiO_2 delivered is 21% which is normal Atmospheric O_2 content

FiO_2 delivery with the Reservoir bag is mathematically modelled, calibrated and verified in the Lab to provide +/- 5% accuracy

Front-panel HMI guides the user in setting the appropriate input O_2 flow rate on the connected O_2 source for a given FiO_2

The mathematical model provides for a possible O_2 concentrator as an O_2 source (purity < 100%)

Online Web-accessible FiO_2 calculator is also provided for exploration purposes



O_2 Flow Rate Calculator
RESPIMATIC 100

Required Incoming O_2 Flow

6.8 (litres/min)

Altitude: 3000 feet ▾

Desired VT(ml) Desired RR(bpm)




O_2 Source Purity(%) Desired FiO_2 (%)



Alarms and Safety Features


- Enforcement of Pre-use checks
- Distinction between Alarms, Warnings and Notifications
- Maintenance Breaths till Alarm situation rectified
 - Max Pressure Alarm
 - Pressure Loss Alarm
 - System Temperature Alarm
 - Mechanical jam Alarm
 - Sensor failure Alarm
 - Breathing Circuit Failure Alarm
 - Detect coughing fits
 - Inconsistent input parameters
 - And many more ...



Power Consumption



An online, sine-wave, external battery UPS recommended to continue operation during power outages



50 AH Car battery is sufficient to run the system for 5+ hours



100 AH Tubular battery is sufficient to run the system for 10+ hours



Input Voltage	180-250 V
Power Consumption	< 100 Watts

Respimatic Testing Process

Timing, Flow and Pressure Checks
for all combinations of various settings



Automated Testing for all combinations of VT, RR, IE, PEEP, PS with randomized spontaneous breath triggers



Automated testing of full day runs checked for timing within 1% of theoretical expectations




All testing so far with a simple test lung
Next step needs a more sophisticated test lung


Future Roadmap



A version using a
mini turbine




A version using
compressed air




Neo-natal &
Invasive Version



Add Vitals (SPO2,
Pulse etc.) to the
Dashboard



Dual Limb
Breathing circuit



More powerful
Microcontroller

About Us

*Evangelizing
Community Service
through Technology*



An extension of the community service initiative that delivered PARC
(Pilani AtmaNirbhar Resource Center*)
www.atmanirbhar.org



Establish pathways that students from various institutes can follow to serve our rural community



Establish a stepping-stone for more sophisticated systems targeting the bottom of the pyramid in future

* The term AtmaNirbhar has been in use by PARC since 2008, and our URL!

Thank you

Backup Details

Market Analysis

Sample Waveforms



The Problem

Scarcity & Affordability

- India has amongst the lowest per capita ICU beds in the World*
- 1.46 beds / 1000 people*
- 3.65 ICU beds / Lakh people*
- Only half of ICU beds are equipped with Ventilators
- A meagre 1.8 Ventilators for one lakh people*
- ICU ventilators are expensive equipment
- Unaffordable in remote clinics
- Ventilator Ambulances are
- Either non-existent except in few major urban centres,
- Or beyond the reach of majority of population

** As per April 2020 Study by Center for Disease Dynamics, Economics & Policy at Princeton University, USA*

The Problem Skilled Practitioners

ICU Ventilators require highly skilled manpower to operate and monitor

Ventilators lie unutilised due to shortage of doctors

TODAY

TIMES NEWS NETWORK

Jaisalmer: The state government has given 17 ventilators, including 12 to Jaisalmer and five to Pokhran government hospitals. However, they were lying unutilised due to lack of doctors. The serious patients are being referred to Jodhpur and a large number of corona patients have died while undergoing treatment in Jodhpur.

On Tuesday, there were 42 fresh cases of Covid-19 in Jaisalmer district. On Sunday, 54 cases were reported. The condition of some patients is serious as they are being referred to Jodhpur which is a five hours journey from Jaisalmer.

tion. The main reason behind referring serious patients to Jodhpur is non-availability of ventilator facility and posts of main doctors are lying vacant.

There are only two physicians in the Jaisalmer hospital of which duty of one of the doctors is to take

व्यवस्थाओं को कोरोना: सरकार ने दिए वेंटीलेटर, आधे से अधिक इंस्टाल नहीं किए

प्रदेश में मरीजों को सासे उखड़ रहीं,
यहाँ स्टोर में 'शो-पीस' बने वेंटीलेटर

कहीं पर्याप्त प्रशिक्षित
स्टाफ ही नहीं
पत्रिका न्यूज़ नेटवर्क
लिकम



दूर के स्थित भरतीय अस्पताल के रुपरे में रखे वेंटीलेटर।

स्वास्थ्य सेवन में वेंटीलेटर में वाई में लाता दिए हैं। तीन को भी एक-दो जो बढ़कर 42 हो गए हैं। इनमें से कुछ फोलिंग वेंटीलेटर हैं। वेंटीलेटर पर मरीज का रखने के लिए आर्सीपी का प्रशिक्षित स्टाफ चाहिए और नियन्त्रण के लिए वेंटीलेटर इंस्टाल कर दिया जाएगा।

अस्पताल में कोरोनावायरस के स्क्रॉप्ट के बाद ऐसे वेंटीलेटर की कमी होने पर सरकार ने वेंटीलेटर उत्पादन की दिए, लेकिन अस्पताल प्रशासन की लापरवाही के चलते जाति चुने में आई वेंटीलेटर अभी तक स्टोर में ही है। बही, चिरोड़ियां वेंटीलेटर के जिला अस्पताल में पर्याप्त प्रशिक्षित स्टाफ नहीं होने से परेशानी आ रही है तथा यहाँ भी 25 वेंटीलेटर इंस्टाल ही नहीं किया गया।

दूसरे स्थित डेढ़जय भरतीय अस्पताल में कोरोनावायरस के स्क्रॉप्ट राहत की वेंटीलेटर ही है। बाद में प्रशासनीय राहत की वेंटीलेटर में 30 नए वेंटीलेटर खेल गए थे। वर्षामास में यहाँ 27 वेंटीलेटर उत्पादन किया जा रहे हैं। 15 इंस्टाल तक नहीं होने से गए।

अ. रमेशन जयपाल, एसोसिएट प्रोफेसर, मैडिकल कॉलेज, दूर

tional doctor to run the ventilator whereas there is need of minimum two to three physicians, cardiologists etc.

Jaisalmer collector Ashish Modi said that all the 17 ventilators in the district are in operational condition and oxygen and other resources are available. He said posts of cardiologist, physician are lying vacant due to which ventilators cannot be used for corona patients. Serious patients are referred to Jodhpur on time and Jodhpur divisional commissioner Dr Samit Sharma is monitoring the situation, he said. Jaisalmer government hospital PMO Dr VK Verma said that ventilator



MUCH TO IMPROVE

23 वेंटीलेटर को इंस्टाल होने का इंतजार

कोरोना भरतीय अस्पताल के बढ़ने के बावजूद और अस्पताल में 52 वेंटीलेटर हैं। और भी बाकी हैं। वहाँ, कोरोना के संबंधित दिए पर्याप्त कामियां हैं।

25 वेंटीलेटर इंस्टाल ही नहीं किए गए

विरोड़ियां वेंटीलेटर में वाई में लाता दिए हैं। तीन को भी एक-दो जो बढ़कर 42 हो गए हैं। इनमें से कुछ फोलिंग वेंटीलेटर हैं। वेंटीलेटर पर मरीज का रखने के लिए आर्सीपी का प्रशिक्षित स्टाफ चाहिए और नियन्त्रण के लिए वेंटीलेटर इंस्टाल कर दिया जाएगा।

दूसरे स्थित विरोड़ियां वेंटीलेटर की विक्रियाकांक्षी नियन्त्रण की विक्रियाकांक्षी वाहियों द्वारा अवृत्ति दी गयी है। वहाँ कुल 42 वेंटीलेटर में से अभी भी पांच पुराने एवं सात नए वेंटीलेटर ही काम में आ रहे हैं। 25 तो इंस्टाल नहीं किये गए।



कुवेरा (नागोर). स्थानीय संसद हमगम कैमियाल की अनुशासन पर शारीर के राजकीय सम्बद्धियक व्यास्त्य केन्द्र को मिले पैटेंटेट वेंटीलेटर की सीएसी के स्टोर में रख दिया गया है। विविधकालिमियों ने बताया कि पैटेंटेट वेंटीलेटर विक्री में रख दिया गया है। तथा जैसे ही जरूरत पहुँची। वार्ड में लेकर काम में ले लिया जाएगा।



Market Need

**25 to 50 ventilators per lakh people
in developed countries**

**1.8 ventilators per lakh people*
(India)**



**Even less for lower income
developing and underdeveloped
countries***

* As per April 2020 Study by Center for Disease Dynamics,
Economics & Policy at Princeton University, USA



Respimatic Sample Waveforms

*Pressure , Flow and
Spontaneous Breaths*



Breath Synchronization in ACV Mode

Tidal volume is delivered at regular intervals T_i .


If spontaneous breath occurs during this interval at time T_s , the set tidal volume is delivered in response, and the next mandatory breath is set for $(T_s + T_i)$.

Again, if another spontaneous breath is detected before $T_s + T_i$, a mandatory breath scheduled T_i time in future.

Breath Syncing in Synchronized AC Mode

There is no sync-window – the next mandatory breath is always rescheduled after a spontaneous breath

Example below: Tidal Volume = 500ml Respiration Rate = 15 bpm



Breath Synchronization in SIMV Mode


Tidal volume is delivered at regular intervals T_i .

If spontaneous breath occurs during this interval at time T_s , a pressure supported breath is delivered in response, and the next mandatory breath is set for $(T_s + 0.9 T_i)$.

This is to allow more spontaneous breaths but to trigger a mandatory VC breath if spontaneous breath is not detected within 90% of the mandatory breath interval.

Breath Syncing in SIMV mode

There is a sync-window – the next mandatory breath is rescheduled only if spontaneous breath within the sync-window



Exceptional Mandatory Breath in PSV Mode


Tidal volume is delivered at regular intervals T_i .

If spontaneous breath occurs during this interval at time T_s , a pressure supported breath is delivered in response.

The next mandatory VC breath is not delayed at all unless it is too close to the spontaneous breath (within 20% of T_i).

This avoids breath stacking issues.

Mandatory Breath (Warning) in PSV mode
There is a quiet-window – the next mandatory breath is delivered only if no spontaneous breath within the quiet-window



Expectations from the session

*Gather feedback on overall
product concept.*

*Do we have the right
price/feature tradeoff ?*

*Is the feature set missing
something critical ?*

*Need access to a sophisticated
test lung for final calibration and
tuning.*

How can you help?

